

Premium Sterile Surgical Blade

KIATO⁺ PLUS



The smart choice for sharp blades

- Carbon Steel Blades made of Swiss steel grade 125Cr1.
- Stainless steel blades made of premium Swedish Steel.
- Individually packed in easy peel color coded foil and boxes.
- Fits all makes of handles. Conforms to ISO 7740.
- Sterilized by Gamma Radiation to a minimum dosage of 2.5 MRads.
- Blades conform to BS 2982 and Medical Device Directive 93/42/EEC.



100 blades per box, 10 box per case, 10 case per master shipping carton.



Size	10	10A	11	11/11P	12	12D	13	15	15/15T	15C
Carbon Steel	18010	18009	18011	18031	18012	18025	18013	18035	18015	18026
Stainless Steel	18110	18109	18111	18131	18112	18125	18113	18135	18115	18126



Size	18	19	20	21	22	22A	23	24	25	26	36
Carbon Steel	18027	18019	18020	18021	18022	18030	18023	18024	18028	18033	18029
Stainless Steel	18127	18119	18120	18121	18122	-	18123	18124	18128	18133	18129

Blade size not to scale.

EC REP
MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany



aditya
DISPOMED

Aditya Dispomed Products Pvt. Ltd.
Plot No. 19 Sector 6 IMT Manesar
Gurgaon Haryana 122050 India
Tel. +91 124 47 64 900
Fax +91 124 47 64 911
info@adityadispomed.com
www.adityadispomed.com

A subsidiary of Kehr Surgical Pvt. Ltd.